

Anti-PROX1 antibody

Cat. No.	ml262082
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-PROX1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human PROX1
Reactivity	Human, Mouse
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	PROX1
Full name	prospero homeobox 1
Synonyms	
Swissprot	Q92786

Target Background

The protein encoded by this gene is a member of the homeobox transcription factor family. Members of this family contain a homeobox domain that consists of a 60-amino acid helix-turn-helix structure that binds DNA and RNA. The protein encoded by this gene is conserved across vertebrates and may play an essential role during development.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 83 kDa

Positive control: 231 cells

Recommended dilution: 200-1000

Gel: 6% SDS-PAGE


Lysate: 40 µg

Lane: 231 cells

Primary antibody: ml262082 (PROX1 Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn