

订购热线: 4008-898-798

Anti-RCC1 antibody

Cat. No. ml262116

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-RCC1 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human RCC1

ReactivityHuman, MouseContent0.3 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol RCC1

Full name regulator of chromosome condensation 1

Synonyms CHC1; RCC1-I; SNHG3-RCC1

Swissprot P18754

Target Background

RCC1 (regulator of chromosome condensation 1), increases the rate at which Ran exchanges GDP for GTP. Ran GAP1 opposes the effects of RCC1 by increasing the rate at which Ran hydrolyzes GTP to GDP. A protein designated Ran BP1 has no intrinsic GAP activity and functions as a GEF inhibitor, deactivating RCC1 and thereby indirectly increasing the ratio of Ran GDP to Ran GTP. Ran BP2 has been proposed as the Ran GTP docking site at the periphery of the nuclear pore complex.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:45 kDa Positive control:Raji cells

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE Lysate: 40 µg Lane: Raji cells

Primary antibody: ml262116(RCC1 Antibody) at dilution

1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds

Da 95 — 72 — 55 — 36 — 28 —

ELISA

Recommended dilution:1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn