

# Anti-TRIM45 antibody

Cat. No.	ml262369
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol
Product overview	
Description	Anti-TRIM45 rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Synthetic peptide of human TRIM45
Reactivity	Human, Mouse
Content	0.4 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
	Antigen affinity purification
Target information	oroc
Symbol	TRIM45
Full name	tripartite motif containing 45
C	200 <sup></sup>
	<i>,</i>



Synonyms RNF99

Swissprot Q9H8W5

## **Target Background**

This gene encodes a member of the tripartite motif family. The encoded protein may function as a transcriptional repressor of the mitogen-activated protein kinase pathway. Alternatively spliced transcript variants have been described.



订购热线: 4008-898-798

#### Applications


### Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human liver cancer

Recommended dilution: 25-100





Good elisa Mil producer:

The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml262369(TRIM45 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### Western blotting

Predicted band size:64 kDa


Positive control:231, Jurkat, A549, hela and 293T cells


Recommended dilution: 500-2000

Predicted cell location: Cytoplasm

Positive control: Human thyroid cancer

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml262369(TRIM45 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)




#### Gel: 6%SDS-PAGE

Exposure time: 1 minute

Lysate: 40 µg

Lane 1-5: 231 cells, Jurkat cells, A549 cells, hela cells, 293T cells Primary antibody: ml262369(TRIM45 Antibody) at dilution 1/400 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution



## ELISA

Recommended dilution: 2000-5000

# 联系电话: 4008-898-798, 021-61725725

# 联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn