

Anti-TGM7 antibody

Cat. No.	ml262406
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-TGM7 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human TGM7
Reactivity	Human, Mouse
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	TGM7
Full name	transglutaminase 7
Synonyms	TGMZ
Swissprot	Q96PF1

Target Background

Transglutaminases (TGM; EC 2.3.2.13) are a family of structurally and functionally related enzymes that stabilize protein assemblies through the formation of gamma-glutamyl-epsilon lysine crosslinks. For additional background information on transglutaminases, see TGM1 (MIM 190195).

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 80 kDa

Positive control: Mouse eyes tissue

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane: Mouse eyes tissue

Primary antibody: ml262406(TGM7 Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn