

Anti-TSPYL1 antibody

Cat. No.	ml262425
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-TSPYL1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human TSPYL1
Reactivity	Human, Mouse
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	TSPYL1
Full name	TSPY-like 1
Synonyms	TSPYL
Swissprot	Q9H0U9

Target Background

The protein encoded by this gene is found in the nucleolus and is similar to that of a family of genes on the Y-chromosome. This gene is intronless. Defects in this gene are a cause of sudden infant death with dysgenesis of the testes syndrome (SIDDT).

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 49 kDa

Positive control: Mouse brain and heart tissue

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-2: Mouse brain tissue,
Mouse heart tissue

Primary antibody:
ml262425(TSPYL1 Antibody) at
dilution 1/200

Secondary antibody: Goat anti
rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn

mlbio 酶联生物
Good elisakit producers
订购热线: 4008-898-798