

Anti-NUP98 antibody

Cat. No.	ml262465
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-NUP98 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human NUP98
Reactivity	Human, Mouse, Rat
Content	1.1 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	NUP98
Full name	nucleoporin 98kDa
Synonyms	ADIR2; NUP96; NUP196
Swissprot	P52948

Target Background

Signal-mediated nuclear import and export proceed through the nuclear pore complex (NPC), which is comprised of approximately 50 unique proteins collectively known as nucleoporins. The 98 kDa nucleoporin is generated through a biogenesis pathway that involves synthesis and proteolytic cleavage of a 186 kDa precursor protein. This cleavage results in the 98 kDa nucleoporin as well as a 96 kDa nucleoporin, both of which are localized to the nucleoplasmic side of the NPC.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 198 kDa

Positive control: K562, 293T, 231 and Hela cells

Recommended dilution: 500-2000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-4: K562 cells, 293T cells, 231 cells, Hela cells

Primary antibody: ml262465(NUP98 Antibody) at dilution 1/550

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn