

Anti-CD79A antibody

Cat. No.	ml262503
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-CD79A rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human CD79A
Reactivity	Human, Mouse
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	CD79A
Full name	CD79a molecule, immunoglobulin-associated alpha
Synonyms	IGA; MB-1
Swissprot	P11912

Target Background

The B lymphocyte antigen receptor is a multimeric complex that includes the antigen-specific component, surface immunoglobulin (Ig). Surface Ig non-covalently associates with two other proteins, Ig-alpha and Ig-beta, which are necessary for expression and function of the B-cell antigen receptor. This gene encodes the Ig-alpha protein of the B-cell antigen component. Alternatively spliced transcript variants encoding different isoforms have been described.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 25 kDa

Positive control: Mouse lung and bladder tissue

Recommended dilution: 500-2000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane 1-2: Mouse lung and bladder tissue

Primary antibody: ml262503(CD79A Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn