

订购热线: 4008-898-798

# Anti-FYN antibody

**Cat. No.** ml262504

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-FYN rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human FYN

Reactivity Human, Mouse, Rat

Content0.5 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol FYN

Full name FYN proto-oncogene, Src family tyrosine kinase

Synonyms SLK; SYN; p59-FYN

Swissprot P06241

# **Target Background**

This gene is a member of the protein-tyrosine kinase oncogene family. It encodes a membrane-associated tyrosine kinase that has been implicated in the control of cell growth. The protein associates with the p85 subunit of phosphatidylinositol 3-kinase and interacts with the fyn-binding protein. Alternatively spliced transcript variants encoding distinct isoforms exist.



订购热线: 4008-898-798

# **Applications**

# Western blotting

Predicted band size:61 kDa Positive control:A549, hela and 293T cell Recommended dilution: 500-2000

Gel: 8%SDS-PAGE

Lysate: 40 µg

Lane 1-3: A549, hela and

293T cell

Primary antibody:


ml262504(FYN Antibody) at

dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000

dilution

Exposure time: 10 seconds



# **ELISA**

elisakii producė Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn



订购热线: 4008-898-798