

Anti-KRT20 antibody

Cat. No.	ml262650
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-KRT20 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human KRT20
Reactivity	Human, Mouse, Rat
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	KRT20
Full name	keratin 20, type I
Synonyms	K20; CD20; CK20; CK-20; KRT21
Swissprot	P35900

Target Background

The protein encoded by this gene is a member of the keratin family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. The type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. This cytokeratin is a major cellular protein of mature enterocytes and goblet cells and is specifically expressed in the gastric and intestinal mucosa. The type I cytokeratin genes are clustered in a region of chromosome 17q12-q21.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 48 kDa

Positive control: Human fetal intestines tissue lysate

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: Human fetal intestines tissue lysate

Primary antibody: ml262650(KRT20 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn