

Anti-SLC4A3 antibody

Cat. No.	ml262842
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-SLC4A3 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human SLC4A3
Reactivity	Human, Mouse, Rat
Content	0.5 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	SLC4A3
Full name	solute carrier family 4 (anion exchanger), member 3
Synonyms	AE3; SLC2C
Swissprot	P48751

Target Background

Plasma membrane anion exchange protein of wide distribution. Mediates at least a part of the Cl⁻/HCO₃⁻ exchange in cardiac myocytes. Both BAE3 and CAE3 forms transport Cl⁻.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 136 kDa

Positive control: A431 and TM4 cell, Human fetal brain tissue

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-3: A431 and TM4 cell, Human fetal brain tissue

Primary antibody: ml262842(SLC4A3 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 90 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn