

订购热线: 4008-898-798

# Anti-UXT antibody

**Cat. No.** ml123519

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-UXT rabbit polyclonal antibody

Applications ELISA, IHC

ImmunogenFull length fusion proteinReactivityHuman, Mouse, Rat

Content0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol UXT

Full name ubiquitously-expressed, prefoldin-like chaperone

Synonyms STAP1; ART-27

Swissprot Q9UBK9

### **Target Background**

The protein encoded by this gene functions as a cofactor that modulates androgen receptor-dependent transcription, and also plays a critical role in tumor necrosis factor-induced apoptosis. Expression of this gene may play a role in tumorigenesis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.



订购热线: 4008-898-798


### **Applications**


## **Immunohistochemistry**

Predicted cell location: Nucleus and Cytoplasm

Positive control: Human prostate cancer

Recommended dilution: 20-100





The image on the left is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using ml123519(UXT Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn