

Anti-C18orf8 antibody

Cat. No.	ml124325
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-C18orf8 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human C18orf8
Reactivity	Human
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	C18orf8
Full name	chromosome 18 open reading frame 8
Synonyms	MIC1; Mic-1; HsT2591
Swissprot	Q96DM3

Target Background

This gene encodes a colon cancer associated protein.

订购热线: 4008-898-798



Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human colorectal cancer

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using ml124325(C18orf8 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: $\times 200$)

Predicted cell location: Cytoplasm

Positive control: Human tonsil

Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml124325(C18orf8 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn