

## Anti-CCDC113 antibody

|                 |                                                         |
|-----------------|---------------------------------------------------------|
| <b>Cat. No.</b> | ml124433                                                |
| <b>Package</b>  | 25 µl/100 µl/200 µl                                     |
| <b>Storage</b>  | -20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol |

### Product overview

|                     |                                         |
|---------------------|-----------------------------------------|
| <b>Description</b>  | Anti-CCDC113 rabbit polyclonal antibody |
| <b>Applications</b> | ELISA, IHC                              |
| <b>Immunogen</b>    | Fusion protein of human CCDC113         |
| <b>Reactivity</b>   | Human, Mouse                            |
| <b>Content</b>      | 0.4 mg/ml                               |
| <b>Host species</b> | Rabbit                                  |
| <b>Ig class</b>     | Immunogen-specific rabbit IgG           |
| <b>Purification</b> | Antigen affinity purification           |

### Target information

|                  |                                   |
|------------------|-----------------------------------|
| <b>Symbol</b>    | CCDC113                           |
| <b>Full name</b> | coiled-coil domain containing 113 |
| <b>Synonyms</b>  | HSPC065                           |
| <b>Swissprot</b> | Q9H0I3                            |

### Target Background

Component of centriolar satellites contributing to primary cilium formation.



## Applications

### Immunohistochemistry

Predicted cell location: Cytoplasm or Nucleus

Positive control: Human ovarian cancer

Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml124433(CCDC113 Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: ×200)

Predicted cell location: Cytoplasm or Nucleus

Positive control: Human esophagus cancer

Recommended dilution: 30-150



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml124433(CCDC113 Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: ×200)

### ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)