

订购热线: 4008-898-798

Anti-CYP20A1 antibody

Cat. No. ml124705

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-CYP20A1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human CYP20A1

ReactivityHumanContent0.54 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol CYP20A1

Full name cytochrome P450 family 20 subfamily A member 1

Synonyms CYP-M Swissprot Q6UW02

Target Background

This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases that catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein lacks one amino acid of the conserved heme binding site. It also lacks the conserved I-helix motif AGX(D,E)T, suggesting that its substrate may carry its own oxygen.




订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm and Cell membrane

Positive control: Human liver cancer Recommended dilution: 20-100





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml124705(CYP20A1 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn