

订购热线: 4008-898-798

Anti-FYTTD1 antibody

Cat. No. ml124858

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-FYTTD1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human FYTTD1

Reactivity Human, Mouse, Rat

Content0.9 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol FYTTD1

Full name forty-two-three domain containing

Synonyms UIF

Swissprot Q96QD9

Target Background

Required for mRNA export from the nucleus to the cytoplasm. Acts as an adapter that uses the DDX39B/UAP56-NFX1 pathway to ensure efficient mRNA export and delivering to the nuclear pore. Associates with spliced and unspliced mRNAs simultaneously with ALYREF/THOC4.



订购热线: 4008-898-798

Applications


Immunohistochemistry

Predicted cell location: Nucleus

Positive control: Human ovarian cancer

Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml124858(FYTTD1 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn