

## Anti-WHRN antibody

<b>Cat. No.</b>	ml263706
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-WHRN rabbit polyclonal antibody
<b>Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic peptide of human WHRN
<b>Reactivity</b>	Human, Mouse
<b>Content</b>	0.66 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	WHRN
<b>Full name</b>	whirlin
<b>Synonyms</b>	WI; CIP98; USH2D; DFNB31; PDZD7B
<b>Swissprot</b>	Q9P202

### Target Background

This gene is thought to function in the organization and stabilization of stereocilia elongation and actin cytoskeletal assembly, based on studies of the related mouse gene. Mutations in this gene have been associated with autosomal recessive non-syndromic deafness and Usher Syndrome. Alternative splicing of this gene results in multiple transcript variants encoding different isoforms.

订购热线: 4008-898-798

#### Applications

##### Western blotting

Predicted band size: 97 kDa

Positive control: 293T, PC-3 and TM4 cell lysates

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40  $\mu$ g

Lane 1-3: 293T, PC-3 and TM4 cell lysates

Primary antibody: ml263706(WHRN Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 minutes



#### ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)