

# Anti-ITSN2 antibody

**Cat. No.** ml264262

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-ITSN2 rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Synthetic peptide of human ITSN2

Reactivity Human, Mouse

Content 0.78 mg/ml

Host species Rabbit

Ig class Immunogen-specific rabbit IgG

Purification Antigen affinity purification

**Target information** 

Symbol ITSN2

Full name intersectin 2



**Synonyms** SWA; SWAP; SH3D1B; SH3P18; PRO2015

Swissprot Q9NZM3

## **Target Background**

This gene encodes a cytoplasmic protein which contains SH3 domains. This protein is a member of a family of proteins involved in clathrin-mediated endocytosis. Intersectin 2 is thought to regulate the formation of clathrin-coated vesicles and also may function in the induction of T cell antigen receptor (TCR) endocytosis.




#### **Applications**

## **Immunohistochemistry**

Predicted cell location: Cytoplasm

Positive control: Human cervical cancer

Recommended dilution: 25-100





Good elisakit producere

The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using ml264262(ITSN2 Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### Western blotting

Predicted band size:193 kDa

Positive control:Hela cell lysate

Recommended dilution: 200-1000



Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane: Hela cell lysate

Primary antibody: ml264262(ITSN2 Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 90 seconds



### **ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn