

订购热线: 4008-898-798

Anti-IRX5 antibody

Cat. No. ml264277

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-IRX5 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human IRX5

ReactivityHuman, MouseContent1.08 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol IRX5

Full nameiroquois homeobox 5SynonymsHMMS; IRXB2; IRX-2a

Swissprot P78411

Target Background

This gene encodes a member of the iroquois homeobox gene family, which are involved in several embryonic developmental processes. Knockout mice lacking this gene show that it is required for retinal cone bipolar cell differentiation, and that it negatively regulates potassium channel gene expression in the heart to ensure coordinated cardiac repolarization. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:50 kDa

Positive control:Mouse lung tissue lysate

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE Lysate: 40 μg

Lane: Mouse lung tissue lysate

Primary antibody: ml264277(IRX5 Antibody) at dilution

1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 3 seconds

kDa 250— 130— 95— 72— 55—

ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn