

兔抗 A4GNT 多克隆抗体

中文名称：兔抗 A4GNT 多克隆抗体

英文名称：Anti-A4GNT rabbit polyclonal antibody

别名：alpha4GnT

抗原：A4GNT

储存：冷冻（-20℃）w

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein from the glycosyltransferase 32 family. The enzyme catalyzes the transfer of N-acetyl glucosamine (GlcNAc) to core 2 branched O-glycans. It forms a unique glycan, GlcNAc α 1 $\rightarrow$ 4Gal β 2 $\rightarrow$ Man α 6 and is largely associated with the Golgi apparatus membrane.
Applications:	ELISA, IHC
Name of antibody:	A4GNT
Immunogen:	Fusion protein of human A4GNT
Full name:	alpha-1,4-N-acetylglucosaminyltransferase
Synonyms :	alpha4GnT

SwissProt:	Q9UNA3
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and human lung cancer
IHC Recommend dilution:	25-100

