

兔抗 ABCA13 多克隆抗体

中文名称: 兔抗 ABCA13 多克隆抗体

英文名称: Anti-ABCA13 rabbit polyclonal antibody

储 存: 冷冻(-20℃) 避光

抗原: ABCA13

宿 主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:

In human, the ATP-binding cassette (ABC) family of trans membrane transporters has at least 48 genes and 7 gene subfamilies. This gene is a member of ABC gene subfamil y A (ABCA). Genes within the ABCA family typically encod e several thousand amino acids. Like other ABC transmem brane transporter proteins, this protein has 12 or more transmembrane alpha-helix domains that likely arrange to f orm a single central chamber with multiple substrate bind ing sites. It is also predicted to have two large extracellul ar domains and two nucleotide binding domains as is typical for ABCA proteins. Alternative splice variants have be en described but their biological validity has not been demonstrated.



Applications:	IF
Name of antibody:	ABCA13
Immunogen:	Synthesized peptide derived from internal of human ABC A13.
Full name:	ATP-binding cassette, sub-family A (ABC1), member 13
Synonyms:	0
SwissProt:	Q86UQ4
IF positive control:	A549 cells
IF Recommend dilution:	100-500

