

兔抗 ABCB5 多克隆抗体

中文名称：兔抗 ABCB5 多克隆抗体

英文名称：Anti-ABCB5 rabbit polyclonal antibody

别名：ATP binding cassette subfamily B member 5; ABCB5beta; EST422562; ABCB5alpha

抗原：ABCB5

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	ABCB5 belongs to the ATP-binding cassette (ABC) transporter superfamily of integral membrane proteins. These proteins participate in ATP-dependent transmembrane transport of structurally diverse molecules ranging from small ions, sugars, and peptides to more complex organic molecules.
Applications:	ELISA, WB
Name of antibody:	ABCB5
Immunogen:	Fusion protein of human ABCB5
Full name:	ATP binding cassette subfamily B member 5
Synonyms:	ABCB5beta; EST422562; ABCB5alpha

SwissProt:	Q2M3G0
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	139 kDa
WB Positive control:	293T cell lysate
WB Recommended dilution:	500-2000

