

兔抗 ABCF3 多克隆抗体

中文名称：兔抗 ABCF3 多克隆抗体

英文名称：Anti-ABCF3 rabbit polyclonal antibody

别 名：EST201864

抗 原：ABCF3

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This protein may display an antiviral effect against flaviviruses such as west Nile virus (WNV) in the presence of OAS1B. ATP-binding cassette (ABC) transporters regulate the transport of a variety of physiologic substrates. But this protein lacks transmembrane domains and is probably not involved in transport.
Applications:	ELISA, IHC
Name of antibody:	ABCF3
Immunogen:	Fusion protein of human ABCF3
Full name:	ATP-binding cassette, sub-family F (GCN20), member 3

Synonyms :	EST201864
SwissProt:	Q9NUQ8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and human lung cancer
IHC Recommend dilution:	25-100

