

## 兔抗 ABHD2 多克隆抗体

中文名称: 兔抗 ABHD2 多克隆抗体

英文名称: Anti-ABHD2 rabbit polyclonal antibody

别名: abhydrolase domain containing 2; HS1-2; LABH2; PHPS1-2

储存: 冷冻 (-20℃)

抗原: ABHD2

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein containing an alpha/beta hydrolase fold, which is a catalytic domain found in a wide range of enzymes. The encoded protein is an acylglycerol lipase that catalyzes the hydrolysis of endocannabinoid arachidonoylglycerol from the cell membrane. This leads to activation of the sperm calcium channel CatSper, which results in sperm activation. Alternative splicing of this gene results in two transcript variants encoding the same protein.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	ABHD2
<b>Immunogen:</b>	Synthetic peptide of human ABHD2
<b>Full name:</b>	abhydrolase domain containing 2
<b>Synonyms:</b>	HS1-2; LABH2; PHPS1-2
<b>SwissProt:</b>	P08910
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer and Human tonsil
<b>IHC Recommend dilution:</b>	40-200



