

兔抗 ABHD6 多克隆抗体

中文名称：兔抗 ABHD6 多克隆抗体

英文名称：Anti-ABHD6 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：ABHD6

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Lipase that preferentially hydrolysis medium-chain saturated monoacylglycerols including 2-arachidonoylglycerol (PubMed:22969151). Through 2-arachidonoylglycerol degradation may regulate endocannabinoid signaling pathways. May also have a lysophosphatidyl lipase activity with a preference for lysophosphatidylglycerol among other lysophospholipids (By similarity).
Applications:	ELISA, IHC
Name of antibody:	ABHD6
Immunogen:	Synthetic peptide of human ABHD6
Full name:	abhydrolase domain containing 6
SwissProt:	Q9BV23
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human esophagus cancer and human tonsil
IHC Recommend dilution:	25-100



