

兔抗 ABHD14B 多克隆抗体

中文名称：兔抗 ABHD14B 多克隆抗体

英文名称：Anti-ABHD14B rabbit polyclonal antibody

别名：CIB; HEL-S-299

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ABHD14B

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Has hydrolase activity towards p-nitrophenyl butyrate (in vitro). May activate transcription.
Applications:	ELISA, WB, IHC
Name of antibody:	ABHD14B
Immunogen:	Fusion protein of human ABHD14B
Full name:	abhydrolase domain containing 14B
Synonyms:	CIB; HEL-S-299
SwissProt:	Q96IU4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human tonsil
IHC Recommend dilution:	50-200

WB Predicted band size:	22 kDa
WB Positive control:	Human fetal liver tissue lysate
WB Recommended dilution:	500-2000



