

## 兔抗 ABL1(Phospho-Thr735) 多克隆抗体

中文名称：兔抗 ABL1(Phospho-Thr735) 多克隆抗体

英文名称：Anti-ABL1(Phospho-Thr735) rabbit polyclonal antibody

别名：ABL; JTK7; p150; c-ABL; v-abl; c-ABL1; bcr/abl

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ABL1(Phospho-Thr735)

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a protooncogene that encodes a protein tyrosine kinase involved in a variety of cellular processes, including cell division, adhesion, differentiation, and response to stress. The activity of the protein is negatively regulated by its SH3 domain, whereby deletion of the region encoding this domain results in an oncogene. The ubiquitously expressed protein has DNA-binding activity that is regulated by CDC2-mediated phosphorylation, suggesting a cell cycle function.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	ABL1(Phospho-Thr735)

<b>Immunogen:</b>	Peptide sequence around phosphorylation site of threonine735 (S-V-T(p)-L-P) derived from Human ABL1.
<b>Full name:</b>	ABL proto-oncogene 1, non-receptor tyrosine kinase
<b>Synonyms :</b>	ABL; JTK7; p150; c-ABL; v-abl; c-ABL1; bcr/abl
<b>SwissProt:</b>	P00519
<b>WB Predicted band size:</b>	123 kDa
<b>WB Positive control:</b>	COS cells lysates treated with EGF
<b>WB Recommended dilution:</b>	500-1000

