

兔抗 ABLIM3 多克隆抗体

中文名称：兔抗 ABLIM3 多克隆抗体

英文名称：Anti-ABLIM3 rabbit polyclonal antibody

别名：actin binding LIM protein family member 3; HMFN1661

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ABLIM3

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the actin-binding LIM (abLIM) family of proteins. These proteins are characterized by an N-terminal LIM domain and a C-terminal dematin-like domain. The encoded protein interacts with actin filaments and may be a component of adherens junctions in several cell types. A variant of this gene may be associated with pain sensitivity in male human patients.
Applications:	ELISA, WB
Name of antibody:	ABLIM3
Immunogen:	Fusion protein of human ABLIM3

Full name:	actin binding LIM protein family member 3
Synonyms:	HMFN1661
SwissProt:	O94929
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	78 kDa
WB Positive control:	PC3 and NIH/3T3 cell lysates
WB Recommended dilution:	500-2000

