

## 兔抗 ACAT2 多克隆抗体

- 中文名称: 兔抗 ACAT2 多克隆抗体
- 英文名称: Anti-ACAT2 rabbit polyclonal antibody
- 相关类别: 一抗
- 储存: 冷冻(-20℃)
- 宿 主: Rabbit
- 抗 原: ACAT2
- 反应种属: Human, Mouse, Rat
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

## 技术规格

Applications: Name of antibody:	n opposite strands of DNA, as well as in opposite tra nscriptional orientation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. ELISA, WB ACAT2
Background:	The product of this gene is an enzyme involved in li pid metabolism, and it encodes cytosolic acetoacetyl- CoA thiolase. This gene shows complementary overla pping with the 3-prime region of the TCP1 gene in both mouse and human. These genes are encoded o



SwissProt:	Q9BWD1
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	41 kDa
WB Positive control:	Jurkat cell lysate
WB Recommended dilution:	500-2000

