

兔抗 ACER3 多克隆抗体

中文名称：兔抗 ACER3 多克隆抗体

英文名称：Anti-ACER3 rabbit polyclonal antibody

别名：APHC; PHCA

抗原：ACER3

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Hydrolyzes only phytoceramide into phytosphingosine and free fatty acid. Does not have reverse activity.
Applications:	WB, IHC, IF
Name of antibody:	ACER3
Immunogen:	Synthesized peptide derived from internal of human P HCA.
Full name:	alkaline ceramidase 3
Synonyms :	APHC; PHCA
SwissProt:	Q9NUN7
IHC positive control:	Human placenta tissue
IHC Recommend dilution:	50-100

WB Predicted band size:	32 kDa
WB Positive control:	COLO cells lysate
WB Recommended dilution:	500-3000
IF positive control:	MCF-7 cells
IF Recommend dilution:	100-500



