

兔抗 ACOT7 多克隆抗体

中文名称：兔抗 ACOT7 多克隆抗体

英文名称：Anti-ACOT7 rabbit polyclonal antibody

别名：acyl-CoA thioesterase 7 ACT; ACH1; BACH; LACH; LACH1; hBACH; CTE-II

抗原：ACOT7

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the acyl coenzyme family. The encoded protein hydrolyzes the CoA thioester of palmitoyl-CoA and other long-chain fatty acids. Decreased expression of this gene may be associated with mesial temporal lobe epilepsy. Alternatively spliced transcript variants encoding distinct isoforms with different subcellular locations have been characterized.
Applications:	ELISA, IHC
Name of antibody:	ACOT7
Immunogen:	Fusion protein of human ACOT7
Full name:	acyl-CoA thioesterase 7

Synonyms:	ACT; ACH1; BACH; LACH; LACH1; hBACH; CTE-II
SwissProt:	O00154
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain
IHC Recommend dilution:	20-100

