

兔抗 ACOX3 多克隆抗体

- 中文名称: 兔抗 ACOX3 多克隆抗体
- 英文名称: Anti-ACOX3 rabbit polyclonal antibody
- 抗 原: ACOX3
- 储存: 冷冻(-20℃)
- 宿 主: Rabbit
- 相关类别: 一抗
- 反应种属: Human
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格



Applications:	ELISA, IHC
Name of antibody:	ACOX3
Immunogen:	Fusion protein of human ACOX3
Full name:	acyl-CoA oxidase 3, pristanoyl
SwissProt:	015254
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and human lung cancer
IHC Recommend dilution:	25-100

