

兔抗 ACP4 多克隆抗体

中文名称：兔抗 ACP4 多克隆抗体

英文名称：Anti-ACP4 rabbit polyclonal antibody

别名：acid phosphatase 4; ACPT; AI1J

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ACP4

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

Acid phosphatases are enzymes capable of hydrolyzing orthophosphoric acid esters in an acid medium. This gene is up-regulated by androgens and is down-regulated by estrogens in the prostate cancer cell line. This gene exhibits a lower level of expression in testicular cancer tissues than in normal tissues. The protein encoded by this gene has structural similarity to prostatic and lysosomal acid phosphatases. Alternatively spliced transcript variants have been described, but their biological validity has not been determined. [provided by RefSeq, Jul 2008]

Applications:	ELISA, IHC
Name of antibody:	ACP4
Immunogen:	Synthetic peptide of human ACP4
Full name:	acid phosphatase 4
Synonyms:	ACPT; AI1J
SwissProt:	Q9BZG2
IHC positive control:	Human thyroid cancer and Human gastric cancer
IHC Recommend dilution:	25-100

