

兔抗 ACP5 多克隆抗体

中文名称：兔抗 ACP5 多克隆抗体

英文名称：Anti-ACP5 rabbit polyclonal antibody

别名：acid phosphatase 5, tartrate resistant; HPAP; TRAP; TRACP5a; TRACP5b; TrATPase

相关类别：一抗

储存：冷冻（-20℃）

抗原：ACP5

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an iron containing glycoprotein which catalyzes the conversion of orthophosphoric monoester to alcohol and orthophosphate. It is the most basic of the acid phosphatases and is the only form not inhibited by L(+)-tartrate.
Applications:	WB, IF
Name of antibody:	ACP5
Immunogen:	Fusion protein of human ACP5
Full name:	acid phosphatase 5, tartrate resistant
Synonyms:	HPAP; TRAP; TRACP5a; TRACP5b; TrATPase
SwissProt:	P13686

WB Predicted band size:	37 kDa
WB Positive control:	BT474, Jurkat cells
WB Recommended dilution:	500-1000
IF positive control:	U2OS cell
IF Recommend dilution:	50-200



