

兔抗 ACTG2 多克隆抗体

中文名称：兔抗 ACTG2 多克隆抗体

英文名称：Anti-ACTG2 rabbit polyclonal antibody

别名：ACT; ACTE; VSCM; ACTA3; ACTL3; ACTSG

储存：冷冻（-20℃）避光

抗原：ACTG2

宿主：Rabbit

相关类别：一抗

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells.
Applications:	WB, IHC
Name of antibody:	ACTG2
Immunogen:	Synthesized peptide derived from human Actin-γ2.
Full name:	actin, gamma 2, smooth muscle, enteric
Synonyms :	ACT; ACTE; VSCM; ACTA3; ACTL3; ACTSG
SwissProt:	P63267
IHC positive control:	Human muscle tissue
IHC Recommend dilution:	50-100

WB Predicted band size:	42 kDa
WB Positive control:	COLO205 cells lysate
WB Recommended dilution:	500-3000



