

兔抗 ACTL8 多克隆抗体

中文名称：兔抗 ACTL8 多克隆抗体

英文名称： Anti-ACTL8 rabbit polyclonal antibody

别 名： CT57

抗 原： ACTL8

储 存： 冷冻（-20℃）

宿 主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This protein belongs to the actin family. Actin-like 8 is a protein in humans that is encoded by the ACTL8 gene. The high levels expression in testis and pancreas and weak expression in placenta. The subcellular location of this protein in cytoplasm, cytoskeleton.
Applications:	ELISA, WB, IHC
Name of antibody:	ACTL8
Immunogen:	Fusion protein of human ACTL8
Full name:	Actin-like 8
Synonyms :	CT57
SwissProt:	Q9H568

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and human esophagus cancer
IHC Recommend dilution:	15-50
WB Predicted band size:	41 kDa
WB Positive control:	Hela and A549 cells
WB Recommended dilution:	500-2000



