

兔抗 ACTN2/3 多克隆抗体

中文名称：兔抗 ACTN2/3 多克隆抗体

英文名称：Anti-ACTN2/3 rabbit polyclonal antibody

别名：CMD1AA

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：ACTN2/3

宿主：Rabbit

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.
Applications:	WB, IHC, IF
Name of antibody:	ACTN2/3
Immunogen:	Synthesized peptide derived from N-terminal of human Actin α -2/3.
Full name:	actinin, alpha 2/3
Synonyms :	CMD1AA
SwissProt:	P35609/Q08043
IHC positive control:	Human skeletal muscle tissue

IHC Recommend dilution:	50-100
WB Predicted band size:	104 kDa
WB Positive control:	Hela cells, HepG2 cells and Huvec cells lysates
WB Recommended dilution:	500-3000
IF positive control:	Hela cells
IF Recommend dilution:	100-500



