

## 兔抗 ACTR3 多克隆抗体

中文名称：兔抗 ACTR3 多克隆抗体

英文名称：Anti-ACTR3 rabbit polyclonal antibody

别名：ARP3

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ACTR3

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The specific function of this gene has not yet been determined; however, the protein it encodes is known to be a major constituent of the ARP2/3 complex. This complex is located at the cell surface and is essential to cell shape and motility through lamellipodial actin assembly and protrusion. Three transcript variants encoding two different isoforms have been found for this gene.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	ACTR3
<b>Immunogen:</b>	Fusion protein of human ACTR3
<b>Full name:</b>	actin related protein 3

<b>Synonyms:</b>	ARP3
<b>SwissProt:</b>	P61158
<b>WB Predicted band size:</b>	47 kDa
<b>WB Positive control:</b>	k562, HL60, Hela, Jurkat, Cos1 and Cos7 cells, mouse spleen and liver tissue
<b>WB Recommended dilution:</b>	500-2000

