

兔抗 ACYP1 多克隆抗体

中文名称：兔抗 ACYP1 多克隆抗体

英文名称：Anti-ACYP1 rabbit polyclonal antibody

别名：ACYPE

相关类别：一抗

抗原：ACYP1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the acylphosphatase family. The encoded protein is a small cytosolic enzyme that catalyzes the hydrolysis of the carboxyl-phosphate bond of acylphosphates. Two isoenzymes have been isolated and described based on their tissue localization: erythrocyte (common) type acylphosphatase encoded by this gene, and muscle type acylphosphatase. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	ACYP1
Immunogen:	Fusion protein of human ACYP1

Full name:	acylphosphatase 1
Synonyms:	ACYPE
SwissProt:	P07311
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300
WB Predicted band size:	11 kDa
WB Positive control:	Human fetal brain tissue lysate
WB Recommended dilution:	500-2000



