

兔抗 ADAM10 多克隆抗体

中文名称：兔抗 ADAM10 多克隆抗体

英文名称：Anti-ADAM10 rabbit polyclonal antibody

别名：kuz; AD10; MADM; CD156c; HsT18717

相关类别：一抗

储存：冷冻（-20℃）

抗原：ADAM10

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an ADAM family member that cleaves many proteins including TNF-alpha and E-cadherin. Members of the ADAM family are cell surface proteins with a unique structure possessing both potential adhesion and protease domains. Sheddase, a generic name for the ADAM metalloproteinase, functions primarily to cleave membrane proteins at the cellular surface. Once cleaved, the sheddases release soluble ectodomains with an altered location and function.
Applications:	ELISA, IHC
Name of antibody:	ADAM10

Immunogen:	Fusion protein of human ADAM10
Full name:	ADAM metallopeptidase domain 10
Synonyms :	kuz; AD10; MADM; CD156c; HsT18717
SwissProt:	O14672
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and human liver cancer
IHC Recommend dilution:	25-100

