

## 兔抗 ADAM12 多克隆抗体

- 中文名称：兔抗 ADAM12 多克隆抗体
- 英文名称：Anti-ADAM12 rabbit polyclonal antibody
- 别名：MCMP; MLTN; CAR10; MLTNA; MCMPMLtna; ADAM12-OT1
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：ADAM12
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of a family of proteins that are structurally related to snake venom disintegrins and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. Expression of this gene has been used as a maternal serum marker for pre-natal development. Alternative splicing results in multiple transcript variants encoding different isoforms. Shorter isoforms are secreted, while longer isoforms are membrane-bound form.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	ADAM12
<b>Immunogen:</b>	Synthetic peptide of human ADAM12
<b>Full name:</b>	ADAM metallopeptidase domain 12
<b>Synonyms :</b>	MCMP; MLTN; CAR10; MLTNA; MCMPMltna; ADAM12-OT 1
<b>SwissProt:</b>	O43184
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human prostate cancer
<b>IHC Recommend dilution:</b>	50-200

