

兔抗 ADAMTS13 多克隆抗体

中文名称：兔抗 ADAMTS13 多克隆抗体

英文名称：Anti-ADAMTS13 rabbit polyclonal antibody

别 名：VWF_{CP}; C9orf8; vWF-CP; ADAM-TS13; ADAMTS-13

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：ADAMTS13

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of a family of proteins containing several distinct regions, including a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. The enzyme encoded by this gene specifically cleaves von Willebrand Factor (vWF). Defects in this gene are associated with thrombotic thrombocytopenic purpura. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	ADAMTS13

Immunogen:	Synthetic peptide of human ADAMTS13
Full name:	ADAM metalloproteinase with thrombospondin type 1 motif, 13
Synonyms :	VWFCP; C9orf8; vWF-CP; ADAM-TS13; ADAMTS-13
SwissProt:	Q76LX8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and human gastric cancer
IHC Recommend dilution:	25-100



