

## 兔抗 ADD1 (Phospho-Ser726) 多克隆抗体

中文名称：兔抗 ADD1 (Phospho-Ser726) 多克隆抗体

英文名称：Anti-ADD1 (Phospho-Ser726) rabbit polyclonal antibody

别名：ADDA

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ADD1 (Phospho-Ser726)

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

Adducins are a family of cytoskeleton proteins encoded by three genes (alpha, beta, gamma). Adducin is a heterodimeric protein that consists of related subunits, which are produced from distinct genes but share a similar structure. Alpha- and beta-adducin include a protease-resistant N-terminal region and a protease-sensitive, hydrophilic C-terminal region. Alpha- and gamma-adducins are ubiquitously expressed. In contrast, beta-adducin is expressed at high levels in brain and hematopoietic tissues. Adducin binds with high affinity to Ca<sup>2+</sup>/calmodulin and is a substrate for protein kinases A and C. Alternative splicing results

	s in multiple variants encoding distinct isoforms; however, not all variants have been fully described.
<b>Applications:</b>	WB, IHC, IF
<b>Name of antibody:</b>	ADD1 (Phospho-Ser726)
<b>Immunogen:</b>	Synthetic peptide of human ADD1 (Phospho-Ser726)
<b>Full name:</b>	adducin 1 (alpha) (Phospho-Ser726)
<b>Synonyms :</b>	ADDA
<b>SwissProt:</b>	P35611
<b>IHC positive control:</b>	Human lung carcinoma
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	130 kDa
<b>WB Positive control:</b>	Hela cells
<b>WB Recommended dilution:</b>	500-1000
<b>IF Positive control:</b>	Hela cells
<b>IF Recommended dilution</b>	100-200



