

兔抗 ADIPOR2 多克隆抗体

中文名称：兔抗 ADIPOR2 多克隆抗体

英文名称：Anti-ADIPOR2 rabbit polyclonal antibody

别名：PAQR2; ACDCR2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ADIPOR2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The adiponectin receptors, ADIPOR1 (MIM 607945) and ADIPOR2, serve as receptors for globular and full-length adiponectin (MIM 605441) and mediate increased AMPK (see MIM 602739) and PPAR-alpha (PPARA; MIM 170998) ligand activities, as well as fatty acid oxidation and glucose uptake by adiponectin .
Applications:	ELISA, IHC
Name of antibody:	ADIPOR2
Immunogen:	Synthetic peptide of human ADIPOR2
Full name:	adiponectin receptor 2
Synonyms :	PAQR2; ACDCR2

SwissProt:	Q86V24
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human lung cancer
IHC Recommend dilution:	25-100

