

## 兔抗 ADPGK 多克隆抗体

中文名称：兔抗 ADPGK 多克隆抗体

英文名称：Anti-ADPGK rabbit polyclonal antibody

别名：ADP-GK; 2610017G09Rik

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：ADPGK

宿主：Rabbit

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

|                          |                                                                                                                                                                                |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>       | Catalyzes the phosphorylation of D-glucose to D-glucose 6-phosphate using ADP as the phosphate donor. GDP and CDP can replace ADP, but with reduced efficiency. By similarity. |
| <b>Applications:</b>     | WB                                                                                                                                                                             |
| <b>Name of antibody:</b> | ADPGK                                                                                                                                                                          |
| <b>Immunogen:</b>        | Synthesized peptide derived from internal of human ADPGK.                                                                                                                      |
| <b>Full name:</b>        | ADP-dependent glucokinase                                                                                                                                                      |
| <b>Synonyms :</b>        | ADP-GK; 2610017G09Rik                                                                                                                                                          |
| <b>SwissProt:</b>        | Q9BRR6                                                                                                                                                                         |

|                                 |                     |
|---------------------------------|---------------------|
| <b>WB Predicted band size:</b>  | 54 kDa              |
| <b>WB Positive control:</b>     | Jurkat cells lysate |
| <b>WB Recommended dilution:</b> | 500-3000            |

