

兔抗 ADRBK1 多克隆抗体

中文名称：兔抗 ADRBK1 多克隆抗体

英文名称：Anti-ADRBK1 rabbit polyclonal antibody

别名：GRK2; BARK1; BETA-ARK1

相关类别：一抗

抗原：ADRBK1

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	The product of this gene phosphorylates the beta-2-adrenergic receptor and appears to mediate agonist-specific desensitization observed at high agonist concentrations. This protein is an ubiquitous cytosolic enzyme that specifically phosphorylates the activated form of the beta-adrenergic and related G-protein-coupled receptors. A abnormal coupling of beta-adrenergic receptor to G protein is involved in the pathogenesis of the failing heart.
Applications:	WB
Name of antibody:	ADRBK1
Immunogen:	Synthesized peptide derived from internal of human AR

	BK1.
Full name:	adrenergic, beta, receptor kinase 1
Synonyms :	GRK2; BARK1; BETA-ARK1
SwissProt:	P25098
WB Predicted band size:	80 kDa
WB Positive control:	Huvec cells lysate
WB Recommended dilution:	500-3000

