

兔抗 CHD3 多克隆抗体

中文名称：兔抗 CHD3 多克隆抗体

英文名称：Anti-CHD3 rabbit polyclonal antibody

级 别：ZFH; Mi-2a; Mi2-ALPHA

抗 原：CHD3

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the CHD family of proteins which are characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. This protein is one of the components of a histone deacetylase complex referred to as the Mi-2/NuRD complex which participates in the remodeling of chromatin by deacetylating histones. Chromatin remodeling is essential for many processes including transcription. Autoantibodies against this protein are found in a subset of patients with dermatomyositis. Three alternatively spliced transcripts encoding different isoforms have been

	described.
Applications:	ELISA, IHC
Name of antibody:	CHD3
Immunogen:	Synthetic peptide of human CHD3
Full name:	chromodomain helicase DNA binding protein 3
Synonyms:	ZFH; Mi-2a; Mi2-ALPHA
SwissProt:	Q12873
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	20-100

