

## 兔抗 CHRNA6 多克隆抗体

中文名称：兔抗 CHRNA6 多克隆抗体

英文名称： Anti-CHRNA6 rabbit polyclonal antibody

别名： CHNRA6

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CHRNA6

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes an alpha subunit of neuronal nicotinic acetylcholine receptors. These receptors consist of five subunits and function as ion channels involved in neurotransmission. The encoded protein is a subunit of neuronal nicotinic acetylcholine receptors that mediate dopaminergic neurotransmission and are activated by acetylcholine and exogenous nicotine. Alternatively spliced transcript variants have been observed for this gene. Single nucleotide polymorphisms in this gene have been associated with both nicotine and alcohol dependence.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	CHRNA6
<b>Immunogen:</b>	Synthetic peptide of human CHRNA6
<b>Full name:</b>	cholinergic receptor nicotinic alpha 6 subunit
<b>Synonyms:</b>	CHNRA6
<b>SwissProt:</b>	Q15825
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer and human ovarian cancer
<b>IHC Recommend dilution:</b>	10-50



