

兔抗 AMOTL2 多克隆抗体

中文名称：兔抗 AMOTL2 多克隆抗体

英文名称：Anti-AMOTL2 rabbit polyclonal antibody

别名：LCCP

相关类别：一抗

抗原：AMOTL2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Angiomotin is a protein that binds angiotensin, a circulating inhibitor of the formation of new blood vessels (angiogenesis). Angiomotin mediates angiotensin inhibition of endothelial cell migration and tube formation in vitro. The protein encoded by this gene is related to angiotensin and is a member of the motilin protein family.
Applications:	ELISA, IHC
Name of antibody:	AMOTL2
Immunogen:	Fusion protein of human AMOTL2
Full name:	angiotensin like 2
Synonyms :	LCCP

SwissProt:	Q9Y2J4
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

