

## 兔抗 CIDEB 多克隆抗体

中文名称：兔抗 CIDEB 多克隆抗体

英文名称：Anti-CIDEB rabbit polyclonal antibody

储 存：冷冻（-20℃）

抗 原：CIDEB

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This protein can activate apoptosis. It is inhibited by DFFB and interacts with DFFA and DFFB. It has high levels of expression in liver and small intestine, and has lower levels of expression in colon, kidney and spleen. CideB specifically interacts with VLDL structural protein, apolipoprotein B100 (apoB100), but not with albumin, a PTG cargo protein. It can interact with COPII components, Sar1 and Sec24.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CIDEB
<b>Immunogen:</b>	Fusion protein of human CIDEB
<b>Full name:</b>	cell death-inducing DFFA-like effector b
<b>SwissProt:</b>	Q9UHD4
<b>ELISA Recommended dilution:</b>	2000-10000

<b>IHC positive control:</b>	Human thyroid cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	50-200

