

兔抗 ANGPTL7 多克隆抗体

中文名称：兔抗 ANGPTL7 多克隆抗体

英文名称：Anti-ANGPTL7 rabbit polyclonal antibody

别名：AngX; CDT6; dJ647M16.1; RP4-647M16.2

相关类别：一抗

抗原：ANGPTL7

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Angiopietin-related protein 7 is a protein that in humans is encoded by the ANGPTL7 gene, a member of the ANGPTL family. Highly and specifically expressed in the cornea where is confined to the stromal layer.
Applications:	ELISA, IHC
Name of antibody:	ANGPTL7
Immunogen:	Fusion protein of human ANGPTL7
Full name:	angiopietin-like 7
Synonyms :	AngX; CDT6; dJ647M16.1; RP4-647M16.2
SwissProt:	O43827

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	25-100

